

Anti-LGALSL antibody

Cat. No.	ml264131
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-LGALSL rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human LGALSL
Reactivity	Human, Mouse
Content	0.78 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	LGALSL
Full name	galectin like
Synonyms	GRP; HSPC159
Swissprot	Q3ZCW2

Target Background

Does not bind lactose, and may not bind carbohydrates.

订购热线: 4008-898-798



Applications

Immunohistochemistry

Predicted cell location: Nucleus and Cytoplasm

Positive control: Human ovarian cancer

Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml264131(LGALSL Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: $\times 200$)

Predicted cell location: Nucleus and Cytoplasm

Positive control: Human esophagus cancer

Recommended dilution: 30-150



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml264131(LGALSL Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn