

Anti-GSC antibody

Cat. No.	ml264135
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-GSC rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human GSC
Reactivity	Human, Mouse
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	GSC
Full name	goosecoid homeobox
Synonyms	SAMS
Swissprot	P56915

Target Background

This gene encodes a member of the bicoid subfamily of the paired (PRD) homeobox family of proteins. The encoded protein acts as a transcription factor and may be autoregulatory. A similar protein in mice plays a role in craniofacial and rib cage development during embryogenesis.

订购热线: 4008-898-798


Applications

Immunohistochemistry


Predicted cell location: Nucleus

Positive control: Human breast cancer

Recommended dilution: 20-100




The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml264135(GSC Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: $\times 200$)




Predicted cell location: Nucleus

Positive control: Human esophagus cancer

Recommended dilution: 20-100



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml264135(GSC Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: $\times 200$)



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn