

订购热线: 4008-898-798

# Anti-IL36B antibody

**Cat. No.** ml264242

Package 25  $\mu$ l/100  $\mu$ l/200  $\mu$ l

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-IL36B rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human IL36B

ReactivityHumanContent1.08 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol IL36B

Full name interleukin 36 beta

Synonyms FIL1; FIL1H; IL1F8; IL1H2; IL-1F8; IL-1H2; IL1-ETA; FIL1-(ETA); FIL1-(ETA)

Swissprot Q9NZH7

### **Target Background**

The protein encoded by this gene is a member of the interleukin 1 cytokine family. Protein structure modeling indicated that this cytokine may contain a 12-stranded beta-trefoil structure that is conserved between IL1A (IL-A alpha) and IL1B (IL-1 beta). This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. Two alternatively spliced transcript variants encoding distinct isoforms have been reported.





订购热线: 4008-898-798

### **Applications**

## **Immunohistochemistry**

Predicted cell location: Secreted Positive control: Human tonsil Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml264242(IL36B Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn