

Anti-UPK3B antibody

Cat. No.	ml120020
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-UPK3B rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Fusion protein of human UPK3B
Reactivity	Human
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	UPK3B
Full name	Uroplakin-3b
Synonyms	FLJ32198, MGC10902, p35, UPIIIb
Swissprot	Q9BT76

Target Background

UPK3B is a minor component of the apical plaques of mammalian urothelium that binds and dimerizes with uroplakin-1b (UPK1B, MIM 602380), one of the major conserved urothelium membrane proteins. The other major conserved integral membrane proteins of urothelial plaques are UPK1A (MIM 611557), UPK2 (MIM 611558), and UPK3A (MIM 611559)

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm or Nucleus

Positive control: Human brain

Recommended dilution: 50-200




The image is immunohistochemistry of paraffin-embedded Human brain tissue using ml120020(UPK3B Antibody) at dilution 1/25. (Original magnification: $\times 200$)

Predicted cell location: Cytoplasm or Nucleus

Positive control: Human tonsil

Recommended dilution: 50-200



The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml120020(UPK3B Antibody) at dilution 1/25. (Original magnification: $\times 200$)

Western blotting

Predicted band size: 34 kDa

Positive control: 293T cells

Recommended dilution: 1000-4000

Gel: 10% SDS-PAGE


Lysate: 40 μ g

Lane: 293T cells

Primary antibody: ml120020(UPK3B Antibody) at dilution 1/1000

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 500-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn