

订购热线: 4008-898-798

Anti-YWHAE antibody

Cat. No. ml120040

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-YWHAE rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Fusion protein of human YWHAE

Reactivity Human, Mouse, Rat

Content0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol YWHAE

Full name Tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, epsilon

polypeptide

Synonyms MDS, MDCR, KCIP-1, 14-3-3E

Swissprot P62258

Target Background

This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 100% identical to the mouse ortholog. It interacts with CDC25 phosphatases, RAF1 and IRS1 proteins, suggesting its role in diverse biochemical activities related to signal transduction, such as cell division and regulation of insulin sensitivity. It has also been implicated in the pathogenesis of small cell lung cancer. Two transcript variants, one protein-coding and the other non-protein-coding, have been found for this gene.





订购热线: 4008-898-798

Applications

Immunohistochemistry

Predicted cell location: Cytoplasm, Nucleus Positive control: Human colon cancer Recommended dilution: 25-150





The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using ml120040(YWHAE Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: ×200)

Western blotting

Predicted band size:29 kDa

Positive control:HT-29, 293T, Hela cells and Mouse brain tissue

Recommended dilution: 500-2000

Gel: 10%SDS-PAGE


Lysate: 30 µg

Lane 1-4: HT29 cells, 293T cells, Hela cells, Mouse brain tissue

Primary antibody: ml120040(YWHAE Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 second



ELISA

Recommended dilution: 500-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio cn@yeah.net

网址: www.mlbio.cn