

订购热线: 4008-898-798

Anti-ANXA2 antibody

Cat. No. ml120065

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-ANXA2 rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Fusion protein of human ANXA2

Reactivity Human, Mouse, Rat

Content0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol ANXA2
Full name Annexin A2

Synonyms P36, ANX2, LIP2, LPC2, CAL1H, LPC2D, ANX2L4, PAP-IV

Swissprot P07355

Target Background

This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. This protein functions as an autocrine factor which heightens osteoclast formation and bone resorption. This gene has three pseudogenes located on chromosomes 4, 9 and 10, respectively. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.





订购热线: 4008-898-798

Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml120065(ANXA2 Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: ×200)

Western blotting

Predicted band size:38 kDa Positive control:HepG2 cells

Recommended dilution: 1000-5000


Gel: 8%SDS-PAGE

Lysate: 30 µg Lane: HepG2 cells

Primary antibody: ml120065(ANXA2 Antibody) at dilution 1/640

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 15 seconds



ELISA

Recommended dilution: 1000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio cn@yeah.net

网址: www.mlbio.cn