

Anti-FHIT antibody

Cat. No.	ml120643
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-FHIT rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Fusion protein of human FHIT
Reactivity	Human, Mouse, Rat
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	FHIT
Full name	fragile histidine triad
Synonyms	FRA3B, AP3Aase
Swissprot	P49789

Target Background

This gene, a member of the histidine triad gene family, encodes a diadenosine 5',5''-P₁,P₃-triphosphate hydrolase involved in purine metabolism. The gene encompasses the common fragile site FRA3B on chromosome 3, where carcinogen-induced damage can lead to translocations and aberrant transcripts of this gene. In fact, aberrant transcripts from this gene have been found in about half of all esophageal, stomach, and colon carcinomas. Alternatively spliced transcript variants have been found for this gene.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human thyroid cancer

Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml120643(FHIT Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: ×200)

Western blotting

Predicted band size: 17 kDa

Positive control: Mouse brain and human fetal kidney tissue

Recommended dilution: 500-2000

Gel: 8+12% SDS-PAGE


Lysate: 30 µg

Lane 1-2: Mouse brain tissue, human fetal kidney tissue

Primary antibody: ml120643(FHIT Antibody) at dilution 1/350

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution: 1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn