

Anti-COX19 antibody

Cat. No.	ml120121
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-COX19 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human COX19
Reactivity	Human, Mouse
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	COX19
Full name	COX19 cytochrome c oxidase assembly factor
Synonyms	MGC104475
Swissprot	Q49B96

Target Background

COX19 encodes a cytochrome c oxidase (COX)-assembly protein. The *S. cerevisiae* Cox19 protein may play a role in metal transport to the mitochondrial intermembrane space and assembly of complex IV of the mitochondrial respiratory chain.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human liver cancer

Recommended dilution: 50-200




The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml120121(COX19 Antibody) at dilution 1/50. (Original magnification: ×200)

Predicted cell location: Cytoplasm

Positive control: Human thyroid cancer

Recommended dilution: 50-200



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml120121(COX19 Antibody) at dilution 1/50. (Original magnification: ×200)

ELISA

Recommended dilution: 1000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn