

Anti-IL13RA2 antibody

Cat. No.	ml120167
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-IL13RA2 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human IL13RA2
Reactivity	Human
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	IL13RA2
Full name	interleukin 13 receptor, alpha 2
Synonyms	CT19; IL-13R; IL13BP; CD213A2
Swissprot	Q14627

Target Background

The protein encoded by this gene is closely related to IL13RA1, a subunit of the interleukin 13 receptor complex. This protein binds IL13 with high affinity, but lacks cytoplasmic domain, and does not appear to function as a signal mediator. It is reported to play a role in the internalization of IL13.

订购热线: 4008-898-798


Applications

Immunohistochemistry


Predicted cell location: Cytoplasm

Positive control: Human lymphoma

Recommended dilution: 15-50




The image on the left is immunohistochemistry of paraffin-embedded Human lymphoma tissue using ml120167(IL13RA2 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: $\times 200$)




Predicted cell location: Cytoplasm

Positive control: Human cervical cancer

Recommended dilution: 15-50



The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using ml120167(IL13RA2 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: $\times 200$)



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn