

订购热线: 4008-898-798

Anti-IL1R2 antibody

Cat. No. ml120173

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-IL1R2 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human IL1R2

ReactivityHumanContent0.5 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol IL1R2

Full name Interleukin-1 receptor type 2

Synonyms IL1RB, CD121b, CDw121b, IL-1R-2, IL-1RT2, IL-1RT-2

Swissprot P27930

Target Background

The protein encoded by this gene is a cytokine receptor that belongs to the interleukin 1 receptor family. This protein binds interleukin alpha (IL1A), interleukin beta (IL1B), and interleukin 1 receptor, type I(IL1R1/IL1RA), and acts as a decoy receptor that inhibits the activity of its ligands. Interleukin 4 (IL4) is reported to antagonize the activity of interleukin 1 by inducing the expression and release of this cytokine. This gene and three other genes form a cytokine receptor gene cluster on chromosome 2q12. Alternative splicing results in multiple transcript variants and protein isoforms. Alternative splicing produces both membrane-bound and soluble proteins. A soluble protein is also produced by proteolytic cleavage.




订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cell membrane, Cytoplasm

Positive control: Human lung cancer Recommended dilution: 25-200





The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml120173(IL1R2 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn