

## Anti-MERTK antibody

 Cat. No.
 ml120440

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-MERTK rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human MERTK
Reactivity	Human, Mouse, Rat
Content	0.4 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	MERTK
Full name	c-mer proto-oncogene tyrosine kinase

,00<sup>c</sup>

c-mer proto-oncogene tyrosine kinase MER, RP38, c-mer Q12866

## **Target Background**

Synonyms

Swissprot


This gene is a member of the MER/AXL/TYRO3 receptor kinase family and encodes a transmembrane protein with two fibronectin type-III domains, two Ig-like C2-type (immunoglobulin-like) domains, and one tyrosine kinase domain. Mutations in this gene have been associated with disruption of the retinal pigment epithelium (RPE) phagocytosis pathway and onset of autosomal recessive retinitis pigmentosa (RP).



订购热线: 4008-898-798

## Applications Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human cervical cancer Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using ml120440(MERTK Antibody) at dilution 1/13, on the right is treated with fusion protein. (Original magnification: ×200)

## ELISA

Recommended dilution: 1000-5000

- 联系电话: 4008-898-798, 021-61725725
- 联系QQ: 2881505695,2881505696
- 邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn