

订购热线: 4008-898-798

Anti-RAB9A antibody

Cat. No. ml122963

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-RAB9A rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Fusion protein of human RAB9A

Reactivity Human, Mouse, Rat

Content0.2 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol RAB9A

Full name RAB9A, member RAS oncogene family

SynonymsRAB9SwissprotP51151

Target Background

Ras-related protein Rab-9A is a protein that in humans is encoded by the RAB9Agene. Rab 9A is a 201 amino acid protein that localizes to the cytoplasmic side of the cell membrane, as well as to the membrane of the Golgi apparatus and the ER, and is involved in the transport of proteins between endosomes and the trans Golgi network.





订购热线: 4008-898-798

Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml122963(RAB9A Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)

Western blotting

Predicted band size:23 kDa

Positive control:Hela, hepg2, RAW264.7 and NIH/3T3 cells

Recommended dilution: 1000-5000


Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane 1-4: Hela cells, hepg2 cells, RAW264.7 cells, NIH/3T3 cells Primary antibody: ml122963(RAB9A Antibody) at dilution 1/250 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution ojisakit pr


Exposure time: 30 seconds

Predicted cell location: Cytoplasm Positive control: Human breast cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using ml122963(RAB9A Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)



ELISA

Recommended dilution: 2000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio cn@yeah.net

网址: www.mlbio.cn