

订购热线: 4008-898-798

Anti-TUBB antibody

Cat. No. ml123070 25 µl/100 µl/200 µl Package Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview Description Applications Immunogen Reactivity Content **Host species** Ig class Purification

Anti-TUBB rabbit polyclonal antibody ELISA, WB, IHC Full length fusion protein Human 0.4 mg/ml Rabbit Immunogen-specific rabbit IgG Antigen affinity purification

Target information	
Symbol	TUBB
Full name	tubulin beta class I
Synonyms	M40; TUBB1; TUBB5; CDCBM6; OK/SW-cl.56
Swissprot	P07437
Target Background	du ^b

Target Background

This gene encodes a beta tubulin protein. This protein forms a dimer with alpha tubulin and acts as a structural component of microtubules. Mutations in this gene cause cortical dysplasia, complex, with other brain malformations 6. Alternative splicing results in multiple splice variants. There are multiple pseudogenes for this gene on chromosomes 1, 6, 7, 8, 9, and Goode 13.




订购热线: 4008-898-798


Applications


Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human cervical cancer Recommended dilution: 25-100



Predicted cell location: Cytoplasm Positive control: Human thyroid cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffinembedded Human cervical cancer tissue using ml123070(TUBB Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml123070(TUBB Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

Western blotting

Predicted band size:50 kDa

Positive control:NIH/3T3, HEPG2, Hela and 293T cell, Mouse brain tissue lysates Recommended dilution: 500-2000

Gel: 8%SDS-PAGE

Lysate: 40 µg

Lane 1-5: NIH/3T3, HEPG2, Hela and 293T cell, Mouse brain tissue

lysates

Primary antibody: ml123070(TUBB Antibody) at dilution 1/300.

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution elisakit P

Exposure time: 10 seconds

kDa 1 2 З 5 130 95 72 55 36-28 17

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio cn@yeah.net

网址: www.mlbio.cn