

订购热线: 4008-898-798

Anti-SLITRK4 antibody

Cat. No. ml123171

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-SLITRK4 rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Fusion protein of human SLITRK4

ReactivityHuman, MouseContent0.1 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol SLITRK4

Full name SLIT and NTRK-like family, member 4

Synonyms

Swissprot Q8IW52

Target Background

This gene encodes a transmembrane protein belonging to the SLITRK family. These family members include two N-terminal leucine-rich repeat domains similar to those found in the axonal growth-controlling protein SLIT, as well as C-terminal regions similar to neurotrophin receptors. Studies of an homologous protein in mouse suggest that this family member functions to suppress neurite outgrowth. Alternative splicing results in multiple transcript variants.




订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human thyroid cancer

Recommended dilution: 10-50





The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml123171(SLITRK4 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: ×200)

Western blotting

Predicted band size:94 kDa Positive control:PC3 cells

Recommended dilution: 200-1000


Gel: 6%SDS-PAGE

Lysate: 40 µg Lane: PC3 cells

don die de la constant de la constan Primary antibody: ml123171(SLITRK4 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 3 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio cn@yeah.net

网址: www.mlbio.cn