

订购热线: 4008-898-798

Anti-GM2A antibody

Cat. No. Package Storage ml123862 25 μl/100 μl/200 μl -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview Description Applications Immunogen Reactivity Content Host species Ig class Purification

Anti-GM2A rabbit polyclonal antibody ELISA, WB, IHC Fusion protein of human GM2A Human, Mouse 0.6 mg/ml Rabbit Immunogen-specific rabbit IgG Antigen affinity purification

Target informationSymbolGM2Full nameGM2SynonymsSAPSwissprotP179

GM2A GM2 ganglioside activator SAP-3; GM2-AP P17900

Target Background

This gene encodes a small glycolipid transport protein which acts as a substrate specific co-factor for the lysosomal enzyme beta-hexosaminidase A. Beta-hexosaminidase A, together with GM2 ganglioside activator, catalyzes the degradation of the ganglioside GM2, and other molecules containing terminal N-acetyl hexosamines. Mutations in this gene result in GM2-gangliosidosis type AB or the AB variant of Tay-Sachs disease. Alternative splicing results in multiple transcript variants.


duce



订购热线: 4008-898-798

Applications

Immunohistochemistry Predicted cell location: Cytoplasm Positive control: Human ovarian cancer Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml123862(GM2A Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: ×200)

Western blotting

0		
redicted band size:21 kDa	netto T	
ositive control:Mouse kidney tissue		
ecommended dilution: 200-1000		
	kDa	
	95-	
Gel: 12%SDS-PAGE	65	
Lysate: 40 µg	36—	
Lane: Mouse kidney tissue	S 28-	
Primary antibody: ml123862(GM2A Antibody) at dilution 1/200)*	
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution	17 —	
Exposure time: 2 seconds		
.x Q ^x	10—	

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

- 邮箱: mlbio_cn@yeah.net
- 网址: www.mlbio.cn