

## Anti-WDFY4 antibody

<b>Cat. No.</b>	ml124089
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-WDFY4 rabbit polyclonal antibody
<b>Applications</b>	ELISA, WB, IHC
<b>Immunogen</b>	Fusion protein of human WDFY4
<b>Reactivity</b>	Human
<b>Content</b>	0.5 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	WDFY4
<b>Full name</b>	WDFY family member 4
<b>Synonyms</b>	C10orf64
<b>Swissprot</b>	Q6ZS81

### Target Background

WDFY4 (WDFY family member 4) is a protein-coding gene. Diseases associated with WDFY4 include chediak-higashi syndrome, and systemic lupus erythematosus.

订购热线: 4008-898-798


#### Applications

##### Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human brain

Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml124089(WDFY4 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification:  $\times 200$ )

##### Western blotting

Predicted band size: 354 kDa

Positive control: Human lung tissue lysate

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40  $\mu$ g

Lane: Human lung tissue lysate

Primary antibody: ml124089(WDFY4 Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 4 minutes



##### ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)