

订购热线: 4008-898-798

# Anti-VASP antibody

**Cat. No.** ml124097

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-VASP rabbit polyclonal antibody

Applications ELISA, WB, IHC

Immunogen Fusion protein of human VASP

ReactivityHuman, MouseContent0.5 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol VASP

Full name vasodilator-stimulated phosphoprotein

**Synonyms** 

Swissprot P50552

#### **Target Background**

Vasodilator-stimulated phosphoprotein (VASP) is a member of the Ena-VASP protein family. Ena-VASP family members contain an EHV1 N-terminal domain that binds proteins containing E/DFPPPXD/E motifs and targets Ena-VASP proteins to focal adhesions. In the mid-region of the protein, family members have a proline-rich domain that binds SH3 and WW domain-containing proteins. Their C-terminal EVH2 domain mediates tetramerization and binds both G and F actin. VASP is associated with filamentous actin formation and likely plays a widespread role in cell adhesion and motility. VASP may also be involved in the intracellular signaling pathways that regulate integrin-extracellular matrix interactions. VASP is regulated by the cyclic nucleotide-dependent kinases PKA and PKG.




订购热线: 4008-898-798

## **Applications**

# **Immunohistochemistry**

Predicted cell location: Cytoplasm Positive control: Human tonsil Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml124097(VASP Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

## Western blotting

Predicted band size:40 kDa

Positive control:A549 and HEPG2 cell lysates

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane 1-2: A549 and HEPG2 cell lysates

Primary antibody: ml124097(VASP Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



### **ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio cn@yeah.net

网址: www.mlbio.cn