

订购热线: 4008-898-798

# Anti-HSPD1 antibody

**Cat. No.** ml120788

Package 25  $\mu$ l/100  $\mu$ l/200  $\mu$ l

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-HSPD1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human HSPD1

Reactivity Human, Mouse, Rat

Content0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol HSPD1

Full name heat shock 60kDa protein 1 (chaperonin)

Synonyms HLD4, CPN60, GROEL, HSP60, HSP65, SPG13, HSP-60, HuCHA60

Swissprot P10809

### **Target Background**

This gene encodes a member of the chaperonin family. The encoded mitochondrial protein may function as a signaling molecule in the innate immune system. This protein is essential for the folding and assembly of newly imported proteins in the mitochondria. This gene is adjacent to a related family member and the region between the 2 genes functions as a bidirectional promoter. Several pseudogenes have been associated with this gene. Two transcript variants encoding the same protein have been identified for this gene. Mutations associated with this gene cause autosomal recessive spastic paraplegia 13.




订购热线: 4008-898-798


### **Applications**

## **Immunohistochemistry**

Predicted cell location: Cytoplasm Positive control: Human ovarian cancer

Recommended dilution: 50-150





The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml120788(HSPD1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn