

Anti-TEKT4 antibody

Cat. No.	ml121332
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-TEKT4 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human TEKT4
Reactivity	Human
Content	0.1 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	TEKT4
Full name	tektin 4
Synonyms	
Swissprot	Q8WW24

Target Background

Structural component of ciliary and flagellar microtubules. Forms filamentous polymers in the walls of ciliary and flagellar microtubules. Found in the abaxial (convex) surface of outer dense fibers in sperm flagella. TEKT4 is necessary for the proper coordinated beating of the sperm flagellum and male reproductive physiology.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human lung cancer

Recommended dilution: 5-20



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml121332(TEKT4 Antibody) at dilution 1/10, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn