

Anti-RARB antibody

Cat. No.	ml120010
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-RARB rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human RARB
Reactivity	Human, Mouse
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	RARB
Full name	retinoic acid receptor, beta
Synonyms	HAP; RRB2; NR1B2
Swissprot	P10826

Target Background

This gene encodes retinoic acid receptor beta, a member of the thyroid-steroid hormone receptor superfamily of nuclear transcriptional regulators. This receptor localizes to the cytoplasm and to subnuclear compartments. It binds retinoic acid, the biologically active form of vitamin A which mediates cellular signalling in embryonic morphogenesis, cell growth and differentiation. It is thought that this protein limits growth of many cell types by regulating gene expression. The gene was first identified in a hepatocellular carcinoma where it flanks a hepatitis B virus integration site. The gene expresses at least two transcript variants; one additional transcript has been described, but its full length nature has not been determined.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:50 kDa

Positive control:Human liver cancer and normal liver tissue

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 50 µg

Lane 1-2: Human liver cancer tissue, Human normal liver tissue

Primary antibody: ml120010(RARB Antibody) at dilution 1/500

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn