

订购热线: 4008-898-798

Anti-ACOT1 antibody

Cat. No. ml121468

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-ACOT1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human ACOT1

ReactivityHumanContent0.3 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol ACOT1

Full nameacyl-CoA thioesterase 1SynonymsACH2; CTE-1; LACH2

Swissprot Q86TX2

Target Background

Acyl-CoA thioesterases such as ACOT1, hydrolyze acyl-CoAs to the free fatty acid and CoA. ACOTs therefore play key roles in maintaining the intracellular ratio between CoA esters of various lipids and free fatty acids. Acyl-CoA thioesterases are a group of enzymes that catalyze the hydrolysis of acyl-CoAs to the free fatty acid and coenzyme A (CoASH), providing the potential to regulate intracellular levels of acyl-CoAs, free fatty acids and CoASH. Active towards fatty acyl-CoA with chain-lengths of C12-C16.



订购热线: 4008-898-798


Applications


Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human esophagus cancer

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml121468(ACOT1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn