

Anti-BSA antibody

Cat. No.	ml120272
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-BSA rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of bovine BSA
Reactivity	Bovine serum albumin
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	BSA
Full name	Bovine serum albumin
Synonyms	
Swissprot	P02769

Target Background

Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca²⁺, Na⁺, K⁺, fatty acids, hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 80% of all plasma zinc.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 69 kDa

Positive control: BSA protein

Recommended dilution: 1000-5000

Gel: 8%SDS-PAGE


Lysate: 0.05 μ g, 0.1 μ g, 0.2 μ g, 0.4 μ g

Lane 1-4: 0.05 μ g, 0.1 μ g, 0.2 μ g, 0.4 μ g BSA

Primary antibody: ml120272(BSA Antibody) at dilution 1/1450

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 second



ELISA

Recommended dilution: 2000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn