

Anti-CCL28 antibody

 Cat. No.
 ml120387

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Anti-CCL28 rabbit polyclonal antibody
ELISA, WB
Fusion protein of human CCL28
Human
0.3 mg/ml
Rabbit
Immunogen-specific rabbit IgG
Antigen affinity purification
CCL28

CCL28 chemokine (C-C motif) ligand 28 MEC, CCK1, SCYA28 Q9NRJ3

Target Background

Full name

Synonyms

Swissprot


This gene belongs to the subfamily of small cytokine CC genes. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterized by two adjacent cysteines. The cytokine encoded by this gene displays chemotactic activity for resting CD4 or CD8 T cells and eosinophils. The product of this gene binds to chemokine receptors CCR3 and CCR10. This chemokine may play a role in the physiology of extracutaneous epithelial tissues, including diverse mucosal organs.



订购热线: 4008-898-798

Applications Western blotting Predicted band size:14 kDa Positive control:Human brain glioma tissue Recommended dilution: 500-2000

Gel: 12%SDS-PAGE Lysate: 40 µg Lane: Human brain glioma tissue Primary antibody: ml120387(CCL28 Antibody) at dilution 1/450 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution Exposure time: 1 minute



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

- 邮箱: mlbio_cn@yeah.net
- 网址: www.mlbio.cn