

订购热线: 4008-898-798

Anti-DKKL1 antibody

Cat. No. ml120555

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-DKKL1 rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Fusion protein of human DKKL1

ReactivityHuman, MouseContent0.6 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol DKKL1

Full name dickkopf like acrosomal protein 1
Synonyms SGY; CT34; SGY-1

Swissprot Q9QZL9

Target Background

The dickkopf protein family interacts with the Wnt signaling pathway and its members are characterized by two conserved cysteine-rich domains. This gene encodes a secreted protein that has low sequence similarity to the dickkopf-3 protein. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:27 kDa

Positive control:Mouse testis tissue lysate

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE

Lysate: 40 µg


Lane: Mouse testis tissue lysate

Primary antibody: ml120555(DKKL1 Antibody) at dilution

1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 second



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn