

Anti-ECE1 antibody

Cat. No.	ml120577
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ECE1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human ECE1
Reactivity	Human, Mouse, Rat
Content	0.9 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ECE1
Full name	endothelin converting enzyme 1
Synonyms	ECE
Swissprot	P42892

Target Background

The protein encoded by this gene is involved in proteolytic processing of endothelin precursors to biologically active peptides. Mutations in this gene are associated with Hirschsprung disease, cardiac defects and autonomic dysfunction. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:87 kDa

Positive control:Mouse heart tissue

Recommended dilution: 1000-5000

Gel: 8%SDS-PAGE


Lysate: 30 μ g

Lane: Mouse heart tissue

Primary antibody: ml120577(ECE1 Antibody) at dilution 1/900

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution:1000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn