

Anti-FOXS1 antibody

Cat. No.	ml120633
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-FOXS1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human FOXS1
Reactivity	Human, Mouse
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	FOXS1
Full name	forkhead box S1
Synonyms	FKHL18, FREAC10
Swissprot	O43638

Target Background

The forkhead family of transcription factors belongs to the winged helix class of DNA-binding proteins. The protein encoded by this intronless gene contains a forkhead domain and is found predominantly in aorta and kidney. The function of the encoded protein is unknown.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 35 kDa

Positive control: Mouse heart tissue

Recommended dilution: 500-2000

Gel: 10+12% SDS-PAGE


Lysate: 40 µg

Lane: Mouse heart tissue

Primary antibody: ml120633 (FOXO1 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn