

Anti-TMPRSS7 antibody

Cat. No.	ml120956
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-TMPRSS7 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human TMPRSS7
Reactivity	Human, Mouse, Rat
Content	0.1 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	TMPRSS7
Full name	transmembrane protease, serine 7
Synonyms	
Swissprot	Q7RTY8

Target Background

Serine protease which preferentially hydrolyzes peptides with Arg at the P1 position. Expressed in brain, ovary, testis, salivary gland, trachea and lung.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 94 kDa

Positive control: A549, hepg2, lovo and A431 cells

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-4: A549 cells, hepg2 cells, lovo cells, A431 cells

Primary antibody: ml120956(TMPRSS7 Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn