

Anti-MLLT3 antibody

Cat. No.	ml121530
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-MLLT3 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human MLLT3
Reactivity	Human, Mouse
Content	0.12 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	MLLT3
Full name	MLLT3, super elongation complex subunit
Synonyms	AF9; YEATS3
Swissprot	P42568

Target Background

Component of the super elongation complex (SEC), a complex required to increase the catalytic rate of RNA polymerase II transcription by suppressing transient pausing by the polymerase at multiple sites along the DNA.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 63 kDa

Positive control: 293T cell lysate

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: 293T cell lysate

Primary antibody: ml121530(MLLT3 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 minutes



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn