

Anti-CRYAA antibody

Cat. No.	ml121589
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-CRYAA rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human CRYAA
Reactivity	Human
Content	0.48 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	CRYAA
Full name	alpha-crystallin A chain
Synonyms	
Swissprot	P02489

Target Background

Contributes to the transparency and refractive index of the lens. Has chaperone-like activity, preventing aggregation of various proteins under a wide range of stress conditions.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 20 kDa

Positive control: Human kidney tissue lysate

Recommended dilution: 200-1000

Gel: 12% SDS-PAGE


Lysate: 40 µg

Lane: Human kidney tissue lysate

Primary antibody: ml121589 (CRYAA Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn