

Anti-MTM1 antibody

Cat. No.	ml121788
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-MTM1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human MTM1
Reactivity	Human, Mouse, Rat
Content	0.42 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	MTM1
Full name	myotubularin 1
Synonyms	CNM; MTMX; XLMTM
Swissprot	Q13496

Target Background

This gene encodes a dual-specificity phosphatase that acts on both phosphotyrosine and phosphoserine. It is required for muscle cell differentiation and mutations in this gene have been identified as being responsible for X-linked myotubular myopathy.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 70 kDa

Positive control: Mouse thymus tissue lysate

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: Mouse thymus tissue lysate

Primary antibody: ml121788(MTM1 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 minutes



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn