

Anti-STK17B antibody

Cat. No.	ml121940
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-STK17B rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human STK17B
Reactivity	Human, Mouse, Rat
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	STK17B
Full name	serine/threonine kinase 17b
Synonyms	DRAK2
Swissprot	O94768

Target Background

Phosphorylates myosin light chains . Acts as a positive regulator of apoptosis. Highly expressed in placenta, lung, pancreas. Lower levels in heart, brain, liver, skeletal muscle and kidney.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 42 kDa

Positive control: A549 cell lysate

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: A549 cell lysate

Primary antibody: ml121940(STK17B Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn