

Anti-CTR9 antibody

Cat. No.	ml161594
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-CTR9 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human CTR9
Reactivity	Human, Mouse
Content	0.96 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	CTR9
Full name	CTR9 homolog, Paf1/RNA polymerase II complex component
Synonyms	TSBP; p150; SH2BP1; p150TSP
Swissprot	Q6PD62

Target Background

The protein encoded by this gene is a component of the PAF1 complex, which associates with RNA polymerase II and functions in transcriptional regulation and elongation. This complex also plays a role in the modification of histones.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 134 kDa

Positive control: K562 cell lysate

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane: K562 cell lysate

Primary antibody: ml161594(CTR9 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10MIN



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn