

订购热线: 4008-898-798

Anti-MB antibody

Cat. No. ml122443

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-MB rabbit polyclonal antibody

Applications ELISA, WB

ImmunogenFull length fusion proteinReactivityHuman, Mouse, Rat

Content 0.3 mg/ml Host species Rabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol MB

Full name myoglobin

Synonyms PVALB; myoglobgin

Swissprot P02144

Target Background

This gene encodes a member of the globin superfamily and is expressed in skeletal and cardiac muscles. The encoded protein is a haemoprotein contributing to intracellular oxygen storage and transcellular facilitated diffusion of oxygen. At least three alternatively spliced transcript variants encoding the same protein have been reported.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:17 kDa

Positive control: Human muscle tissue, Rat heart tissue, Mouse heart tissue and Human heart tissue lysates

Recommended dilution: 500-2000

Gel: 12%SDS-PAGE

Lysate: 40 µg

Lane 1-4: Human muscle tissue, Rat heart tissue, Mouse


heart tissue and Human heart tissue lysates

Primary antibody: ml122443(MB Antibody) at dilution

1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 second



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn