

Anti-GTPBP10 antibody

Cat. No.	ml122454
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-GTPBP10 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human GTPBP10
Reactivity	Human
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	GTPBP10
Full name	GTP-binding protein 10
Synonyms	ObgH2; UG0751c10
Swissprot	A4D1E9

Target Background

Small G proteins, such as GTPBP10, act as molecular switches that play crucial roles in the regulation of fundamental cellular processes such as protein synthesis, nuclear transport, membrane trafficking, and signal transduction

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 43 kDa

Positive control: A549 cells

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: A549 cells

Primary antibody: ml122454(GTPBP10 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn