

Anti-HSD17B7 antibody

Cat. No.	ml122503
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-HSD17B7 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human HSD17B7
Reactivity	Human, Mouse, Rat
Content	0.3 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	HSD17B7
Full name	hydroxysteroid (17-beta) dehydrogenase 7
Synonyms	PRAP; SDR37C1
Swissprot	P56937

Target Background

HSD17B7 encodes an enzyme that functions both as a 17-beta-hydroxysteroid dehydrogenase in the biosynthesis of sex steroids and as a 3-ketosteroid reductase in the biosynthesis of cholesterol.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 38 kDa

Positive control: Human testis tissue

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: Human testis tissue

Primary antibody: ml122503(HSD17B7 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn