

Anti-TXLNA antibody

| | |
|-----------------|---|
| Cat. No. | ml122538 |
| Package | 25 µl/100 µl/200 µl |
| Storage | -20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol |

Product overview

| | |
|---------------------|---------------------------------------|
| Description | Anti-TXLNA rabbit polyclonal antibody |
| Applications | ELISA, WB |
| Immunogen | Fusion protein of human TXLNA |
| Reactivity | Human, Mouse |
| Content | 0.5 mg/ml |
| Host species | Rabbit |
| Ig class | Immunogen-specific rabbit IgG |
| Purification | Antigen affinity purification |

Target information

| | |
|------------------|---------------|
| Symbol | TXLNA |
| Full name | taxilin alpha |
| Synonyms | IL14; TXLN |
| Swissprot | P40222 |

Target Background

α-taxilin is a novel binding partner of the syntaxin family which is implicated in intracellular vesicle trafficking. Through its C-terminal coiled-coil region, α-taxilin interacts with the nascent polypeptide-associated complex (NAC), which acts as a transcriptional coactivator.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 62 kDa

Positive control: A549 and 293T cells

Recommended dilution: 500-2000

Gel: 6% SDS-PAGE


Lysate: 40 µg

Lane 1-2: A549 cells, 293T cells

Primary antibody: ml122538(TXLNA Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn