

Anti-KLHL1 antibody

Cat. No.	ml122589
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-KLHL1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human KLHL1
Reactivity	Human, Mouse
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	KLHL1
Full name	kelch-like family member 1
Synonyms	MRP2
Swissprot	Q9NR64

Target Background

The KLHL1 protein belongs to a family of actin-organizing proteins related to Drosophila Kelch (Nemes et al., 2000 [PubMed 10888605]).

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 83 kDa

Positive control: 293T, PC3 and Jurkat cells

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-3: 293T cells, PC3 cells, Jurkat cells

Primary antibody: ml122589(KLHL1 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn