

Anti-OVAL antibody

Cat. No.	ml122670
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-OVAL rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	OVAL protein
Reactivity	Chicken
Content	0.5 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	OVAL
Full name	ovalbumin (SERPINB14)
Synonyms	SERPINB14
Swissprot	P01012

Target Background

Ovalbumin (abbreviated OVA) is the main protein found in egg white, making up 60-65% of the total protein. Ovalbumin displays sequence and three-dimensional homology to the serpin superfamily, but unlike most serpins it is not a serine protease inhibitor. It is presumed to be a storage protein.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 43 kDa

Positive control: OVAL protein

Recommended dilution: 500-2000

Gel: 10% SDS-PAGE


Lysate: 0.4 μ g, 0.8 μ g

Lane 1-2: 0.4 μ g, 0.8 μ g OVAL protein

Primary antibody: ml122670(OVAL Antibody) at dilution 1/500

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 3 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn