

Anti-NOX5 antibody

Cat. No.	ml122775
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-NOX5 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human NOX5
Reactivity	Human
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	NOX5
Full name	NADPH oxidase, EF-hand calcium binding domain 5
Synonyms	
Swissprot	Q96PH1

Target Background

This gene is predominantly expressed in the testis and lymphocyte-rich areas of spleen and lymph nodes. It encodes a calcium-dependen NADPH oxidase that generates superoxide, and functions as a calcium-dependent proton channel that may regulate redox-dependent processes in lymphocytes and spermatozoa. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 86 kDa

Positive control: Raji, TM4 and PC3 cells

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-3: Raji cells, TM4 cells, PC3 cells

Primary antibody: ml122775(NOX5 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn