

Anti-OTUB2 antibody

Cat. No.	ml122814
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-OTUB2 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human OTUB2
Reactivity	Human, Mouse
Content	0.5 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	OTUB2
Full name	OTU deubiquitinase, ubiquitin aldehyde binding 2
Synonyms	OTB2; OTU2; C14orf137
Swissprot	Q96DC9

Target Background

This gene encodes one of several deubiquitylating enzymes. Ubiquitin modification of proteins is needed for their stability and function; to reverse the process, deubiquitylating enzymes remove ubiquitin. This protein contains an OTU domain and binds Uba1 (ubiquitin aldehyde); an active cysteine protease site is present in the OTU domain.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 27 kDa

Positive control: Human liver cancer tissue

Recommended dilution: 500-2000

Gel: 10% SDS-PAGE


Lysate: 40 µg

Lane: Human liver cancer tissue

Primary antibody: ml122814(OTUB2 Antibody) at dilution 1/800

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn