

Anti-HOXB2 antibody

Cat. No.	ml161752
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-HOXB2 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human HOXB2
Reactivity	Human
Content	0.5 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	HOXB2
Full name	homeobox B2
Synonyms	K8; HOX2; HOX2H; Hox-2.8
Swissprot	P14652

Target Background

This gene is a member of the Antp homeobox family and encodes a nuclear protein with a homeobox DNA-binding domain. It is included in a cluster of homeobox B genes located on chromosome 17. The encoded protein functions as a sequence-specific transcription factor that is involved in development. Increased expression of this gene is associated with pancreatic cancer.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:38 kDa

Positive control:Human bladder cancer tissue

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: Human bladder cancer tissue

Primary antibody: ml161752(HOXB2 Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 minutes



ELISA

Recommended dilution:1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn