

Anti-NUP85 antibody

Cat. No.	ml122847
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-NUP85 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human NUP85
Reactivity	Human, Mouse, Rat
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	NUP85
Full name	nucleoporin 85kDa
Synonyms	Nup75
Swissprot	Q9BW27

Target Background

Bidirectional transport of macromolecules between the cytoplasm and nucleus occurs through nuclear pore complexes (NPCs) embedded in the nuclear envelope. NPCs are composed of subcomplexes, and NUP85 is part of one such subcomplex, Nup107-160.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 75 kDa

Positive control: Hela cells

Recommended dilution: 500-2000

Gel: 6% SDS-PAGE


Lysate: 40 µg

Lane: Hela cells

Primary antibody: ml122847(NUP85 Antibody) at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn