

Anti-S100A11 antibody

Cat. No.	ml123077
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-S100A11 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human S100A11
Reactivity	Human, Mouse, Rat
Content	0.5 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	S100A11
Full name	S100 calcium binding protein A11
Synonyms	MLN70; S100C; HEL-S-43
Swissprot	P31949

Target Background

The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein may function in motility, invasion, and tubulin polymerization.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 12 kDa

Positive control: PC3 cells

Recommended dilution: 1000-5000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane: PC3 cells

Primary antibody: ml123077(S100A11 Antibody) at dilution 1/1000

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



ELISA

Recommended dilution: 2000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn