

订购热线: 4008-898-798

Anti-THRA antibody

Cat. No. ml123324

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-THRA rabbit polyclonal antibody

Applications ELISA, WB

Immunogen Fusion protein of human THRA

Reactivity Human, Mouse, Rat

Content0.4 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol THRA

Full name thyroid hormone receptor, alpha

Synonyms AR7; EAR7; ERBA; CHNG6; ERBA1; NR1A1; THRA1; THRA2; ERB-T-1; c-ERBA-1

Swissprot P10827

Target Background

The protein encoded by this gene is a nuclear hormone receptor for triiodothyronine. It is one of the several receptors for thyroid hormone, and has been shown to mediate the biological activities of thyroid hormone. Knockout studies in mice suggest that the different receptors, while having certain extent of redundancy, may mediate different functions of thyroid hormone.



订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:55 kDa

Positive control:Mouse muscle tissue Recommended dilution: 500-2000

Gel: 8%SDS-PAGE Lysate: 40 μg


Lane: Mouse muscle tissue

Primary antibody: ml123324(THRA Antibody) at dilution

1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution:2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn