

Anti-TAGAP antibody

Cat. No.	ml124808
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-TAGAP rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Full length fusion protein
Reactivity	Human, Mouse
Content	0.78 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	TAGAP
Full name	T-cell activation RhoGTPase activating protein
Synonyms	FKSG15; IDDM21; TAGAP1; ARHGAP47
Swissprot	Q8N103

Target Background

This gene encodes a member of the Rho GTPase-activator protein superfamily. The encoded protein may function as a Rho GTPase-activating protein. Alterations in this gene may be associated with several diseases, including rheumatoid arthritis, celiac disease, and multiple sclerosis. Alternate splicing results in multiple transcript variants encoding distinct isoforms.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 81 kDa

Positive control: A549, Jurkat and HEPG2 cell lysates

Recommended dilution: 500-2000

Gel: 6% SDS-PAGE


Lysate: 40 µg

Lane 1-3: A549, Jurkat and HEPG2 cell lysates

Primary antibody: ml124808 (TAGAP Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn