

Anti-CIZ1 antibody

Cat. No.	ml162097
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-CIZ1 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human CIZ1
Reactivity	Human
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	CIZ1
Full name	CDKN1A interacting zinc finger protein 1
Synonyms	NP94; LSFR1; ZNF356
Swissprot	Q9ULV3

Target Background

The protein encoded by this gene is a zinc finger DNA binding protein that interacts with CIP1, part of a complex with cyclin E. The encoded protein may regulate the cellular localization of CIP1. Several transcript variants encoding different isoforms have been found for this gene.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 100 kDa

Positive control: A549 cells

Recommended dilution: 200-1000

Gel: 6% SDS-PAGE


Lysate: 40 µg

Lane: A549 cells

Primary antibody: ml162097(CIZ1 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn