

Anti-MAGOHB antibody

Cat. No.	ml125248
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-MAGOHB rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Fusion protein of human MAGOHB
Reactivity	Human, Mouse
Content	1.02 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	MAGOHB
Full name	mago homolog B, exon junction complex core component
Synonyms	MGN2; mago; magoh
Swissprot	Q96A72

Target Background

Involved in mRNA splicing and in the nonsense-mediated decay (NMD) pathway.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Nucleus and Cytoplasm

Positive control: Human esophagus cancer

Recommended dilution: 50-300




The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml125248(MAGOHB Antibody) at dilution 1/60. (Original magnification: ×200)

Predicted cell location: Nucleus and Cytoplasm

Positive control: Human thyroid cancer

Recommended dilution: 50-300



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml125248(MAGOHB Antibody) at dilution 1/60. (Original magnification: ×200)

Western blotting

Predicted band size: 17 kDa

Positive control: HL-60, 293T, Hela and Jurkat cell lysates

Recommended dilution: 500-2000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane 1-4: HL-60, 293T, Hela and Jurkat cell lysates

Primary antibody: ml125248(MAGOHB Antibody) at dilution 1/400

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 15 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn