

Anti-C8orf44 antibody

Cat. No.	ml125259
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-C8orf44 rabbit polyclonal antibody
Applications	ELISA, WB, IHC
Immunogen	Fusion protein of human C8orf44
Reactivity	Human
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	C8orf44
Full name	chromosome 8 open reading frame 44
Synonyms	
Swissprot	Q96CB5

Target Background

The function of this protein remains unknown.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human thyroid cancer

Recommended dilution: 40-200




The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml125259(C8orf44 Antibody) at dilution 1/35. (Original magnification: ×200)

Predicted cell location: Cytoplasm

Positive control: Human gastric cancer

Recommended dilution: 40-200



The image is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using ml125259(C8orf44 Antibody) at dilution 1/35. (Original magnification: ×200)

Western blotting

Predicted band size: 18 kDa

Positive control: K562 cell, Human kidney tissue lysates

Recommended dilution: 200-1000

Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane 1-2: K562 cell, Human kidney tissue lysates

Primary antibody: ml125259(C8orf44 Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn