

Anti-RNF5 antibody

Cat. No.	ml125282
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-RNF5 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human RNF5
Reactivity	Human, Mouse, Rat
Content	0.6 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	RNF5
Full name	ring finger protein 5
Synonyms	RMA1; RING5
Swissprot	Q99942

Target Background

The protein encoded by this gene contains a RING finger, which is a motif known to be involved in protein-protein interactions. This protein is a membrane-bound ubiquitin ligase. It can regulate cell motility by targeting paxillin ubiquitination and altering the distribution and localization of paxillin in cytoplasm and cell focal adhesions.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 20 kDa

Positive control: 293T and Jurkat cell lysates

Recommended dilution: 500-2000

Gel: 12% SDS-PAGE


Lysate: 40 µg

Lane 1-2: 293T and Jurkat cell lysates

Primary antibody: ml125282(RNF5 Antibody) at dilution 1/250

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn