

Anti-SCGN antibody

Cat. No.	ml125390
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-SCGN rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Fusion protein of human SCGN
Reactivity	Human, Mouse, Rat
Content	1.8 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	SCGN
Full name	secretagogin, EF-hand calcium binding protein
Synonyms	SEGN; CALBL; SECRET; setagin; DJ501N12.8
Swissprot	O76038

Target Background

The encoded protein is a secreted calcium-binding protein which is found in the cytoplasm. It is related to calbindin D-28K and calretinin. This protein is thought to be involved in KCL-stimulated calcium flux and cell proliferation.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size: 32 kDa

Positive control: Human fetal brain tissue lysate

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 µg

Lane: Human fetal brain tissue lysate

Primary antibody: ml125390(SCGN Antibody) at dilution 1/600

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn